



# AUTOPSY SAW



## **Anatomy Electric Hand Bone Saw**

### **Model No.: SAH**

**Application:** Medicolegal expertizing instrument, surgical instrument  
The autopsy saw can be worked with both AC and DC power.

No damage to cartilage tissue, no splashing during cutting, and it has very strong cutting force on hard substance with plasticity, it cuts very fast.

The saw blades are made of high quality stainless steel, very durable. Simple and convenient in operation, safe and lightweight.

Motor with low temperature, and it doesn't burn.

### **Host parameters:**

Nominal voltage: DC 12V

Unload current:  $\geq 4.5A$

Input power:  $\geq 120W$

Axis swing times: 16000/min

Main unit weight: 1.1kg

Swing angle:  $4^\circ$

### **Control box parameters:**

Input voltage: 100V-240V

Output voltage: 10v-12v

Battery: 86.4Wh

Charging time: 4 hours