

CarnisMini-6KC/6KS/6KT and CarnisMini-6K/7K/10K Mini Centrifuge

Features

1. Composite shell, shatter-resistant
 2. The electronic braking system for rapid sample processing
 3. Quick release rotor tub, allows for easy rotor exchange
 4. CarnisMini-6KC / CarnisMini-6KS / CarnisMini-6KT
 - Two speed options: 4,000 & 6,000rpm
 - . CarnisMini-6KC include 2 quick release interchangeable rotors
 - 8x1.5/2.0ml tube rotor, --4 place PCR Strip rotor
 - . CarnisMini-6KS: 3-in-1 rotor design, high capacity and time-saving
 - 6x1.5/2.0ml tubes, --6x0.5ml tubes, --2-8x0.2ml PCR strip tubes
 - . CarnisMini-6KT: Patented rotor design for 48-well 50ul PCR plate
- Optional rotor: 4x5ml Eppendorf tube



CarnisMini-6KC



CarnisMini-6KS

Patented 3-in-1 rotor design



CarnisMini-6KT




05
CarnisMini-6KT 48-well
50ul PCR plate rotor



01
CarnisMini-6KC quick
release rotor



02
CarnisMini-6KS
3-in-1 rotor



03
CarnisMini-6KT
4x5ml tube rotor



04
CarnisMini-6KC/CarnisMini-6KS
/CarnisMini-6KT
Two speed options



CarnisMini-6K / CarnisMini-7K / CarnisMini-10K

Specification

Model	CarnisMini-6K	CarnisMini-7K	CamisMini-10K	CarnisMini-6KC	CarnisMini-6KS	CamisMini-6KT
Rotor	6x1.5ml/2.0ml rotor 2x8x0.2ml strip rotor	6x1.5ml/2.0ml rotor 2x8x0.2ml strip rotor	6x1.5ml/2.0ml rotor	8x1.5ml/2.0ml rotor 4x8x0.2ml strip rotor	3-in-1 rotor (*)	48-well 50ul PCR rotor or 5x5ml tube rotor
Speed	6,000rpm	7,000rpm	10,000rpm	4,000 & 6,000rpm	4,000 & 6,000rpm	4,000 & 6,000rpm
Max.RCF	2,000 x g	2,300 x g	5,000 x g	1,200 & 2,000 x g	1,200 & 2,000 x g	1,200 & 2,000 x g
Dimensions	175x148x118mm	175x148x118mm	175x148x118mm	135x147x120mm	135x127x120mm	135x147x120mm
Weight	0.7kg	0.7kg	0.7kg	0.7kg	0.7kg	0.7kg
Power	AC220-240V	AC220-240V	AC220-240V	AC220-240V	AC220-240V	AC220-240V

*3-in-rotor ---- 6x1.5/2.0ml, 6x0.5ml tubes and 2-8x0.2ml PCR strip rotor

Ordering information

No.	Code	Description
01	AS-08060-00	CamisMini-6K, AC220-240V, 60Hz
02	AS-08070-00	CamisMini-7K, AC220-240V, 60Hz
03	AS-08080-00	CamisMini-10K, AC220-240V, 60Hz
04	AS-08020-00	CamisMini-6KC, AC220-240V, 60Hz
05	AS-08030-00	CamisMini-6KS, AC220-240V, 60Hz
06	AS-08100-00	CamisMini-6KT, AC220-240V, 60Hz

Accessories

CarnisMini-6K/CarnisMini-7K/CarnisMini-10K		CarnisMini-6KC		CarnisMini-6KS	
AS-08061-01	Allen wrench	AS-08021-01	8 x 1.5/2.0ml rotor	AS-08031-01	6x1.5/2.0ml, 6x0.5ml tubes & 2-8x0.2ml PCR Strip rotor
AS-08061-02	6 x 1.5/2.0ml rotor	AS-08021-02	4 x 8 x 0.2ml strip rotor		
AS-08061-03	2 x 8 x 0.2ml strip rotor	AS-08061-05	Adapters for 0.5/0.6ml tubes	CarnisMini-6KT	
AS-08061-04	Slide rotor	AS-08061-06	Adapters for 0.2ml tubes	AS-08101-01	48 well 50ul PCR plate rotor
AS-08061-05	Adapters for 0.5/0.6ml tubes			AS-08101-02	4x5ml tubes rotor
AS-08061-06	Adapters for 0.2ml tubes			AS-08101-03	Adapter for 5ml normal tubes

THE
DRAKE
BUSINESS GROUP
CORPORATION



Drake Marketing and
Equipment Corporation



Drake Calibration
Laboratory Philippines



Erpel Konstrukt
Konzern Corporation



Urossorsa Marketing and
Equipment Corporation

(02) 8463 4943 | (02) 8463 5120
(+63) 939 911 9277

1052 Benigno Aquino Avenue, Sitio 2
Bahay Pare, Candaba, Pampanga
Philippines, 2013

176 Benigno Aquino Avenue, Sitio 2
Bahay Pare, Candaba, Pampanga
Philippines, 2013

information@drakephilippines.com