

CarnisMPC-25/CarnisMini-P25 Mini-Plate Centrifuge

Features for CarnisMPC-P25

1. Unique design for the horizontal rotor centrifuge (Patented)
2. Short-spin and timing of two operations make the tests more convenient
3. Braking time from max.speed: less than 30sec
4. Timer can be set up to 10 minutes
5. Streamline design with more compact shape and transparent PC lid
6. It can accept below formats:
 - PCR standard plates with skirted, non-skirted
 - Elisa Plates
 - 8x0.2ml PCR strip tubes
 - 0.2ml PCR tubes

CarnisMPC-P25



01

CarnisMPC-P25 Operation



■ PCR plates



■ Elisa plates



■ 2pcs PCR plates with skirt
or without skirt



■ 8 × 0.2ml PCR strip tubes

Specification

Model	CarnisMPC-P25
Max. rotational speed	2,200rpm
Max. RCF	480 x g
Capacity	2pcs PCR plates / Elisa plate / 8x0.2ml PCR strip tubes / 0.2ml PCR tubes
Weight	3.6kg
Dimensions (W x D x H)mm	290 x 360 x 140mm
Power supply	AC220-240V, 60Hz

02

CarnisMPC-P25

Adaptor for PCR plates with skirt and without skirt



03

CamisMPC-P25

Adaptor for 8x0.2ml strip or 0.2ml PCR tubes



Ordering information

No.	Code	Description
01	AS-08090-00	CarnisMPC-P25 Mini-Plate Centrifuge, AC220-240V, 60Hz
02	AS-08091-01	Adaptor for PCR plates with skirt and without skirt
03	AS-08091-02	Adaptor for 8x0.2ml strip or 0.2ml PCR tubes

CarnisMini-P25



Features for CarnisMini-P25

1. Quickly spin down droplets
2. Use before and after thermal cycling to increase PCR yield
3. Accepts skirted, non-skirted and all standard PCR plates
4. Less than 1/4 the size of most plate centrifuges

Specification

Model	CarnisMini-P25
Max. rotational speed	2,500rpm
Max. RCF	500 x g
Capacity	2pcs PCR plate
Weight	2.1kg
Dimensions (W x D x H)mm	190 x 220 x 186
Power supply	AC220-240V, 60Hz

Ordering information

No.	Code	Description
01	AS-08010-00	CarnisMini-P25 Mini-Plate Centrifuge

THE
DRAKE
 BUSINESS GROUP
 CORPORATION



Drake Marketing and Equipment Corporation



Drake Calibration Laboratory Philippines
1052 Benigno Aquino Avenue, Sitio 2 Bahay Pare, Candaba, Pampanga Philippines, 2013



Erpel Konstrukt Konzern Corporation
176 Benigno Aquino Avenue, Sitio 2 Bahay Pare, Candaba, Pampanga Philippines, 2013



Urossorsa Marketing and Equipment Corporation

(02) 8463 4943 | (02) 8463 5120
(+63) 939 911 9277

information@drakephilippines.com