



MICROTOME

KANARIEN Rotary Microtome



Features:

1. It has a clipper-built case, the user can put blade, wax block, or other consumable one the top of case. It is easy.
2. It does not need to lubricate and maintain the microtome frequently.
3. The handwheel on the right can be locked in any position, and the handle can be locked in highest position, it is double safe for changing blade or tissue block.
4. The specimen Orienting system: axes X and Y 8°, axis 180°.
5. It adopts the special blade holder, it is more safe and convenient to use the disposable blade directly. It also can to longer the user life of disposable blade.

Technical Specifications

1. Slice thickness range: 0.5- 60um

2. Setting values : From 0-2 um in 0.5 um –increments

Form 2-10um in 1 um-increments

Form 10-20 um in 2 um-increments

Form 20-60um in 5 um-increments

3. Horizontal specimen stroke :28 mm

4. Vertical specimen stroke : 60mm

5. Slip precision: $\pm 5\%$

6. Maximal slice section: 50 × 45mm

Accessories:

Standard Setting List

- | | |
|----------------------------|-------|
| 1. Disposable Blade Holder | 1 pc |
| 2. Specimen Cassette | 1 pc |
| 3. Hexagonal Wrench | 1 set |
| 4. Brush | 1 pc |

Optional Accessories

Paraffin Wax Holder

Feather Microtome Blade

Package Info:

Dimensions: 620 ×525×455mm,

Weight:45 kgs