

## Pathology Paraffin Block Storage Cabinet

### GEIER-CC01



One piece of paraffin block storage cabinet



Base



GEIER-CC01

Item No.	Details	L x W x H(mm)	Color	Weight(kg)	Unit
	One piece of paraffin block storage cabinet	525*500*380	Ivory	26	Each
	Base	525*500*70	Ivory	1.5	Each
GEIER-CC01	One set of paraffin block storage cabinet	525*500*1590	Ivory	105.5	Each

### FEATURES:

1. Four pieces and a base combine one set of paraffin block storage cabinet
2. Outer dimension: 515\*480\*1510 mm
3. Capacity: each of the 4 piece cabinet holds approximately 4,125 pcs wax block with a total capacity of 16,500 pcs wax block per set.
4. It has the same dimension with our pathology drying slide storage cabinet, pathology slide storage cabinet. So they could be stacked up, saving much space and looked elegant.
5. Material: cold-rolled steel sheet with 0.8 mm thickness.
6. Surface treatment: automatic electrostatic powder coating.
7. Inner structure: each piece has 5 drawers, each drawer has 12 grooves in it, and the groove width is adjustable.
8. The drawer could be pulled out smoothly and steadily, it has automatic homing feature, which could prevent the drawers from falling down accidentally.
9. There's a card slot location for each drawer, so the contents of each drawer could be easily identified.
10. With high elastic shock absorption mat, reduce collisions with cabinet put oneself in another's position to close the drawer, reduce the noise to the minimum.
11. On the top of each cabinet has securing snaps. So they could stack up one by one, saving space and more steadier.
12. Packing volume: 0.4 cbm

## Pathology Paraffin Block Storage Cabinet

### GEIER-CC02



Item No.	Details	L x W x H(mm)	Color	Weight(kg)	Unit
	One piece of paraffin block storage cabinet	515*480*360	Ivory	26	Each
	Base	515*480*70	Ivory	1.5	Each
GEIER-CC02	One set of paraffin block storage cabinet	515*480*1510	Ivory	79.5	Each

### Features:

1. Specifications (mm); 515\*480\*1510 m m
2. Cabinet made using high quality cold-rolled steel plate 1.0m.
3. Drawer: wax block special pumping.
4. Handle: Plastic dark handle.
5. Surface: degreasing, surface conditioning, zinc phosphating, passivation, powder spraying.
6. Slot: special slot for metal mold opening.
7. Standard groove: integrated stamping forming.
8. The overall structure of the combined type, each can store about 12000 pieces of standard wax.
9. Use the standard embedding box drawer storage design.