

Pathology Slide Storage Cabinet

GEIER-SC01



One piece of
slide storage



Base



GEIER-SC03

Item No.	Commodity	L x W x H(mm)	Color	Weight(kg)	Unit
	One piece of slide storage cabinet	515*480*360	Ivory	24	Each
	Base	515*480*70	Ivory	1.5	Each
GEIER-SC03	One set of slide storage cabinet	515*480*1510	Ivory	97.5	Each

FEATURE:

1. Four pieces and a base combine one set of slide storage cabinet
2. Outer dimension: 515*480*1510 mm
3. Capacity: each of the 4 piece cabinet holds approximately 13,750 pcs slides with a total capacity of 55,000 pcs slides per set.
4. It has the same dimension with our pathology drying slide storage cabinet, pathology paraffin wax storage cabinet. So they could be stacked up, saving much space and looked elegant.
5. Material: cold-rolled steel sheet with 0.8 mm thickness.
6. Surface treatment: automatic electrostatic powder coating.
7. Inner structure: each piece has 9 drawers, each drawer is divided into small grooves, and the groove width is adjustable.
8. The drawer could be pulled out smoothly and steadily, it has automatic homing feature, which could prevent the drawers from falling down accidentally.
9. There's a card slot location for each drawer, so the contents of each drawer could be easily identified.
10. With high elastic shock absorption mat, reduce collisions with cabinet put oneself in another's position to close the drawer, reduce the noise to the minimum.
11. On the top of each cabinet has securing snaps. So they could stack up one by one, saving space and more steadier.
12. Packing volume: 0.4 cbm

Pathology Slide Storage Cabinet

GEIER-SC02



One piece of slide storage



GEIER-SC02

Item No.	Commodity	L x W x H(mm)	Color	Weight(kg)	Unit
	One piece of slide storage cabinet	460*550*430	Ivory	23.79	Each
	Base	515*480*70	Ivory	1.5	Each
GEIER-SC02	One set of slide storage cabinet	460*550*1630	Ivory	97	Each

Features:

1. Specification (mm): 450*478*1630mm
2. 72 pumping 12 section, each section 6 pumping, plus a base, each drawer can store the standard glass 900 sections, each section can be stored 6 million or so, cabinet with domestic 1.0m cold rolled steel plate; the drawer slide: ABS special oil drawer slide, built in lock, with anti slip function. Handle with ABS chrome drawer slide by ABS special separator for archiving.
3. The surface of the cabinet: degreasing, surface conditioning, zinc phosphating, passivation, powder spraying.
4. Slot: special slot for metal mold opening.
5. Standard groove: integrated stamping forming.