

WACHTEL-1145 Tissue Flotation Workstation



It is composed in three areas, water Flattening Bath, Slide Drying Bench and a multipurpose Tray.

Features:

1. Water Bath-Slide Dryer is composed in three areas, water Flattening Bath,
2. Slide Drying Bench and a multipurpose Tray that can be used for drying and warming
3. slides. All functions are settable through the control panel: temperature of Water Bath,
4. Drying Bench and Thermal Tray.

Technical Specifications:

1. Dimension: (W-D-H) 550×340×220 mm
2. Dimension of Water Bath: (W-D-H) 135×185×80 mm
3. Temperature:
4. Water Bath: ambient-- 99°C
5. Drying Bench: ambient-- 99°C