



Efficient Diagnostic

HD60EXP

Color Doppler Ultrasonic Diagnostic System



Efficient Diagnosis

THE DRAKE BUSINESS GROUP CORPORATION


 Drake Marketing and Equipment Corporation	 Drake Calibration Laboratory Philippines	 Erpel Konstrukt Konzern Corporation	 Urossorsa Marketing and Equipment Corporation
--	---	--	---

1052 Benigno Aquino Avenue, Sitio 2 Bahay Pare, Candaba, Pampanga Philippines, 2013
176 Benigno Aquino Avenue, Sitio 2 Bahay Pare, Candaba, Pampanga Philippines, 2013
information@drakephilippines.com



(02) 8463 4943 | (02) 8463 5120
(+63) 939 911 9277

HD60EXP system platform

HD60EXP's platform is adopted GPU processing architecture, through multi-line scanning, multi-line receiving, lossless data uploading and floating point operation, the accuracy and speed of signal processing are greatly improved, the image has better whole field uniformity and spatial resolution.





HiScom Spatial Composite Imaging




- Real-time image comparison on one screen
- Support for the color mode
- Support for convex and linear array probes
- Support for harmonic imaging

HiPure Advanced Speckle Noise Suppression

Speckle noise reduction suppression technology suppresses speckles and improves contrast resolution without blurring structural details.



Excellent Imaging in Abdominal, Obstetric, and Gynecological Examination



Intelligent Health Examination Software Hi-Step

Helps physicians finish ultrasound examination more easily



Increasing patients
Faster ultrasound examination required



Insufficient resources and high workload
Higher efficiency of clinical work required



Need for a standard operation procedure
No ultrasound examination item missed



Simple User Operation

Solo ingenuity of ergonomic design, electric lifting, rotating adjustable control panel, omni-directional adjustable screen joint arm, give you a comfortable and intimate working environment.

