



Vortex Mixer

LIFE SCIENCE

Alpinum Prazision Austria GmbH & Co. KG

Description

Mini vortex is the ideal vortex mixer for every lab with its compact size, easy of use, and its economic price point. The 2,500rpm fixed speed ensures there is no sacrificing any mixing power. It requires little working space and is easy-to-use.

Features

1. Touch activated vortex mixer, stable compact design
2. 2,500rpm fixed vortexing speed, powerful vortexing for tubes up to 15ml
3. CarnisMTV-1 is with a power adapter and a battery case



Specification

Model	CarnisMTV-1
Working style	Continues / Short-Spin
Shaking movement	Orbital
Mixing Orbit	3mm
Speed Range	800~2500rpm
Dimensions (mm)	74 x 91 x 84mm
Input Power	DC6V
Weight	0.35Kg

Ordering information

No.	Code	Description
01	AS-02070-00	CarnisMTV-1, DC6V, 0.6W
02	AS-02071-01	DC6V Power Adapter
03	AS-02071-02	Battery box for 4pcs AA batteries

THE
DRAKE
BUSINESS GROUP
C O R P O R A T I O N



Drake Marketing and
Equipment Corporation



Drake Calibration
Laboratory Philippines



Erpel Konstrukt
Konzern Corporation



Urossorsa Marketing and
Equipment Corporation

(02) 8463 4943 | (02) 8463 5120
(+63) 939 911 9277

1052 Benigno Aquino Avenue, Sitio 2
Bahay Pare, Candaba, Pampanga
Philippines, 2013

176 Benigno Aquino Avenue, Sitio 2
Bahay Pare, Candaba, Pampanga
Philippines, 2013

information@drakephilippines.com